

# Aesculap® Power Systems

Catalogue 2011




Aesculap Power Systems

**B|BRAUN**  
SHARING EXPERTISE

# Aesculap® Power Systems

Aesculap – pioneer in the field of surgical power systems.

Aesculap introduced its first surgical electric motor as early as 1904. Further important milestones followed in subsequent decades with the introduction of compressed-air systems in the 1960s, the further development of flexible cable electric motors in the 1970s, the introduction of battery-driven motors in the 1980s and the development of pneumatic and electric high speed systems in the 1990s.



With the current systems Acculan 3Ti, HiLAN XS and microspeed uni, Aesculap offers you customized solutions for virtually all fields of application. Our system range is completed by the innovative Eccos storage system which enables timely, mechanical reprocessing of your Aesculap power systems.

Aesculap Power Systems

....get the speed you need.

## How the catalogue works

### 1. Select motor system

indicates handpiece compatibility

**GD675** **GB586R**  
microspeed uni XS Hi  
High speed motor for all  
Hi-Line XS handpieces  
115 W; 10.000 – 80.000 1/min;  
; 2,6 Ncm; 76 x 18 mm;  
80 g; Lateral

••• ➤ **Hi-Line XS**

consisting parts:  
GB600820 Sterilite spray adaptor Hi  
GB600870 Sterilite spray adaptor XS

Electric

### 2. Select handpiece & Eccos holder

~ 40 mm

**GB756R** **GB614R**  
Hi-Line XS handpiece  
1 ring code;  
spray nozzle GB761SU  
(not included)

article code Eccos  
fixation

Hi-Line XS

### 3. Select Tool Use Burrs & Blades catalogue 017599

**Aesculap® Burrs & Blades**  
Catalogue 2010

**2011**

**BIBRAUN** Medical Systems

Chapter  
Hi-Line XS  
burrs

Hi-Line XS

## Signs and symbols

	<b>Latex-free</b> all products listed in this catalogue are latex-free
	Article code Eccos fixation
	Forward / reverse drive
	Forward drive
	Oscillation frequency
	Oscillaton angle
	Stroke
	Cannulation
	Cable length
	Weight
	Sterile
	Single use
	Packaging Unit
	Required height of tray
	Suitable handpiece-line

# Content

## Power Systems

- **Electric Power Systems**
  - microspeed uni
  - microspeed arthro
  - Acculan 3Ti Electro
  - Elan EC
- **Pneumatic Power Systems**
  - HiLAN XS
- **Battery Power Systems**
  - Acculan 3Ti

## Handpieces

- Hi-Line XS
- micro-Line
- mini-Line

## Accessories

## Processing

- Sterilit Power Systems
- Eccos

## Index



# Power Systems

## Electric Power Systems

- **microspeed uni**
- **microspeed arthro**
- **Acculan 3Ti Electro**
- **Elan EC**

## Pneumatic Power Systems

- **HiLAN XS**

## Battery Power Systems

- **Acculan 3Ti**

# Electric Power Systems

## microspeed uni



### Electric power system for neurosurgery, neurootology and spine

- > Versatile system for small bone surgery and arthroscopy
- > Reliable and durable due to the application of top-grade material and state-of-the art technology
- > Easy operation over touch screen and function control button
- > Compatible with handpiece lines Hi-Line XS, Hi-Line, micro-Line and mini-Line

**GD670****microspeed uni control unit**

with irrigation pump; touch screen;  
power cord not included;  
pump flow 0 - 80 ml/min  
100-240 V AC, 50/60Hz;  
305 x 305 x 175 mm; **8.300 g**  
consisting parts:  
GD412804 bottle holder

**GD671****microspeed uni foot control**

double pedal; protection type IPx8;  
anaesthetic proof class M  
**5.0 m;**  
250 x 225 x 75 mm; **2.390 g**

**GD673****GB689R****microspeed uni motor cable**

with hand control; for use with  
motors GD674, GD676, GD678,  
GD685

**3.0 m; 380 g**

**GD668****microspeed uni foot control**

single pedal; protection type IPx8;  
anaesthetic proof class M  
**5.0 m;**  
250 x 225 x 75 mm; **2.420 g**

**GD672****GB689R****microspeed uni motor cable**

for use with motors GD674, GD676,  
GD678, GD685  
**3.0 m; 360 g**

**GD676****GB578R****microspeed uni Hi 150**

high speed motor for all Hi-Line and  
Hi-Line XS handpieces  
150 W; 10.000 - 80.000 1/min;  
**3.5 Ncm; 99 x 27 mm; 155 g**

**•••► Hi-Line, Hi-Line XS**

consisting parts:  
GB600820 Sterilite spray adaptor Hi

# microspeed uni



**GD675** GB586R

**microspeed uni XS Hi**

high speed motor for all Hi-Line XS handpieces

115 W; 10.000 – 80.000 1/min; ↗; 2.6 Ncm; 76 x 18 mm; 80 g

•••► **Hi-Line XS**

consisting parts:

GB600870 Sterilite spray adaptor XS

GB600820 Sterilite spray adaptor Hi



**GD685**

GB587R

**microspeed uni perforator driver**

to use with cranial perforators with Hudson shank

150 W; 800 – 1.200 1/min; ↗; 2.0 Nm; 155 x 41 mm; 530 g

consisting parts:

GB600820 Sterilite spray adaptor Hi



**GD678** GB578R

**microspeed uni micro 150**

low speed motor for all micro-Line handpieces

150 W; 1.000 – 30.000 1/min; ↗; 7.0 Ncm; 113 x 27 mm; 161 g

•••► **micro-Line**

consisting parts:

GB600810 Sterilite spray adaptor

micro-Line

GB600820 Sterilite spray adaptor Hi



**GD674**

GB578R

**microspeed uni mini 100**

low speed motor for all mini-Line handpieces except of GB385R, GB388R, GB397R, GB398R

100 W; 1.000 – 16.000 1/min; ↗; 7.0 Ncm; 95 x 27 mm; 176 g

•••► **mini-Line**

consisting parts:

GB600820 Sterilite spray adaptor Hi

GB600830 Sterilite spray adaptor

mini-Line



**GD684** GB579R

**microspeed uni pistol handpiece**

low speed pistol for all mini-Line handpieces except of GB373R, GB374R, GB385R, GB397R; oscillation and tapping mode

100 W; 1.000 – 16.000 1/min; ↗; 7.0 Ncm; ↗ 2.5 mm; ↘ 3.0 m; 95 x 27 x 120 mm; 480 g

•••► **mini-Line**

consisting parts:

GB600830 Sterilite spray adaptor

mini-Line

GD684225 K-wire protection sleeve



**GD681**

GB580R

**microspeed uni shaver handpiece**

with hand control; for Aesculap shaver tips; variable suction; automatic tip recognition; oscillation speed adjustable 0.3 – 2.0 seconds

100 W; 400 – 8.000 1/min; ↗; ↘ 3.0 m; 150 x 27 x 37 mm; 510 g

consisting parts:

GB600820 Sterilite spray adaptor Hi

**GD682****GB580R****microspeed uni shaver handpiece**

for Aesculap shaver tips; variable suction; automatic tip recognition; oscillation speed adjustable 0.3 – 2.0 seconds; to be used in combination with foot control GD671

100 W; 400 – 8.000 1/min; ↗;

↙ 3.0 m; 150 x 27 x 37 mm;

510 g

consisting parts:

GB600820 Sterilite spray adaptor Hi

Power cords, irrigation sets and  
trolleys see 'Accessories'

# Electric Power Systems

## microspeed arthro



### Electric power system for arthroscopy and small bone surgery

- > Reliability and durability due to the application of top-grade material and state-of-the-art manufacturing technology
- > Shaver handpieces available with ergonomic two button hand control or with foot control
- > Compact and lightweight construction, with automatic tip recognition
- > Individual adjustment of oscillation frequency and acceleration and stopping characteristics
- > Large range shaver tips and burrs for excellent cutting and ablation performance

**GD665****microspeed arthro control unit**  
touch screen; power cord not included100-240 V AC, 50/60Hz;  
305 x 305 x 175 mm; 5.300 g**GD671****microspeed uni foot control**  
double pedal; protection type IPx8;  
anaesthetic proof class M  
5.0 m;  
250 x 225 x 75 mm; 2.390 g**GD672****GB689R****microspeed uni motor cable**  
for use with motors GD674, GD676,  
GD678, GD685

3.0 m; 360 g

**GD673****GB689R****microspeed uni motor cable**  
with hand control; for use with  
motors GD674, GD676, GD678,  
GD685

3.0 m; 380 g

**GD681****GB580R****microspeed uni shaver handpiece**  
with hand control; for Aesculap  
shaver tips; variable suction;  
automatic tip recognition;  
oscillation speed adjustable 0.3 - 2.0  
seconds100 W; 400 - 8.000 1/min; 3.0 m;  
150 x 27 x 37 mm;  
510 gconsisting parts:  
GB600820 Sterilite spray adaptor Hi**GD682****GB580R****microspeed uni shaver handpiece**  
for Aesculap shaver tips; variable  
suction; automatic tip recognition;  
oscillation speed adjustable 0.3 - 2.0  
seconds; to be used in combination  
with foot control GD671100 W; 400 - 8.000 1/min; 3.0 m;  
150 x 27 x 37 mm;  
510 gconsisting parts:  
GB600820 Sterilite spray adaptor Hi

# microspeed arthro



**GD684** **GB579R**

## microspeed uni pistol handpiece

low speed pistol for all mini-Line handpieces except of GB373R, GB374R, GB385R, GB397R; oscillation and tapping mode  
100 W; 1.000 - 16.000 1/min; ↗; 7.0 Ncm; 2.5 mm; 3.0 m; 95 x 27 x 120 mm; 480 g

### •••► mini-Line

consisting parts:  
GB600830 Sterilite spray adaptor  
mini-Line  
GD684225 K-wire protection sleeve



**GD674** **GB578R**

## microspeed uni mini 100

low speed motor for all mini-Line handpieces except of GB385R, GB388R, GB397R, GB398R  
100 W; 1.000 - 16.000 1/min; ↗; 7.0 Ncm; 95 x 27 mm; 176 g

### •••► mini-Line

consisting parts:  
GB600820 Sterilite spray adaptor Hi  
GB600830 Sterilite spray adaptor  
mini-Line



**GD678** **GB578R**

## microspeed uni micro 150

low speed motor for all micro-Line handpieces  
150 W; 1.000 - 30.000 1/min; ↗; 7.0 Ncm; 113 x 27 mm; 161 g

### •••► micro-Line

consisting parts:  
GB600810 Sterilite spray adaptor  
micro-Line  
GB600820 Sterilite spray adaptor Hi

Power cords, irrigation sets and  
trolleys see 'Accessories'



# Electric Power Systems

## Acculan 3Ti Electro



**Electric power system for orthopaedics, traumatology, cardiothoracic surgery, and dermatology**

- > Powerful - up to 250 W
- > Reliable
- > Easy handling - Titanium housing allows mechanical alkaline reprocessing



**GA871** GB495R  
**Acculan 3Ti Electro small drill**  
with oscillation and tapping mode  
250 W; ↘ 4.0 mm;  
107 x 47 x 180 mm; 1.100 g  
consisting parts:  
GA876 adapter  
GA875 lid  
GA678 sterile funnel  
GA672211 rinsing adaptor  
GA672207 K-wire protection sleeve



**GA873** GB495R  
**Acculan 3Ti Electro oscillating saw**  
250 W; <>>> 0 - 15.000 1/min;  
↙ 4°47'; 107 x 47 x 200 mm;  
1.600 g  
consisting parts:  
GA876 adapter  
GA875 lid  
GA678 sterile funnel



**GA874** GB495R  
**Acculan 3Ti Electro reciprocating saw**  
250 W; <>>> 0 - 15.000 1/min;  
↙ 3.2 mm; 180 x 47 x 200 mm;  
1.540 g  
consisting parts:  
GB440215 cleaning plate  
GA876 adapter  
GA875 lid  
GA678 sterile funnel



**GB436R** GB686R  
**Protection guide**  
for use with GA674, GA874 and saw blade GC761R



**GA870** GB498R  
**Acculan 3Ti Electro dermatome**  
width keyless adjustable; thickness adjustable (1/10 mm steps);  
cutting thickness: 0.2 - 1.2 mm;  
cutting width: max. 78 mm;  
suitable blades GB228R (to be ordered separately)  
250 W; <>>> 0 - 6.500 1/min;  
↙ 3.1 mm; 290 x 110 x 90 mm;  
1.700 g  
consisting parts:  
GA876 adapter  
GA875 lid  
GA678 sterile funnel



**GA877**  
**Acculan 3Ti Electro console**  
power cord not included  
100-240 V AC, 50/60Hz;  
255 x 217 x 90 mm; 2.930 g

# Acculan 3Ti Electro



**GA878** GB486R  
**Acculan 3Ti Electro connection cable**  
for use with GA870, GA871, GA873,  
GA874  
3.0 m; 410 g



**GA875** GB485R  
**Acculan 3Ti Electro lid**



**GA876**  
**Acculan 3Ti Electro adapter**  
120 x 49 x 42 mm; 165 g



**GA678** GB488R  
**Acculan 3Ti sterile funnel**



# Electric Power Systems

## Elan EC



### Universal power system for all disciplines

- > Unlimited power
- > Legendary reliability and durability
- > Easy handling
- > More than 15.000 systems installed worldwide

**GA830****Elan EC control unit**

power cord not included  
100–240 V AC, 50/60Hz;  
305 x 305 x 125 mm; **7.300 g**

**GA835****Elan EC control unit**

with irrigation pump; power cord not included;  
pump flow 0 – 80 ml/min  
100–240 V AC, 50/60Hz;  
305 x 305 x 125 mm; **7.900 g**  
consisting parts:  
GD412804 bottle holder

**GA188****Foot control**

waterproof (IP67) acc. to IEC601  
and DIN40050; anaesthetic proof  
class M  
**5.0 m;**  
130 x 200 x 60 mm; **1.900 g**

**GA176****GB680R****Flexible cable micro**

for all micro-Line handpieces  
0 – 20.000 1/min; **1.8 m;**  
**600 g**

•••► **micro-Line**

**GA173****GB680R****Flexible cable micro**

for all micro-Line handpieces  
0 – 20.000 1/min; **2.3 m;**  
**700 g**

•••► **micro-Line**

Power cords, irrigation sets and  
trolleys see 'Accessories'

# Pneumatic Power Systems

## HiLAN XS



### **Power, reliability and flexibility**

- > HiLAN XS, HiLAN, microLAN - Integration of high and low speed
- > The pneumatic system for all applications in microsurgery

recommended pressure for all pneumatic AESCULAP Systems is 6 - 10 bar, 85 - 145 psi (technical information at 6 bar)

## HiLAN XS



**GA740R** GB610R

### HiLAN XS

pneumatic high speed motor for all Hi-Line XS handpieces  
50 W; 0 - 90.000 1/min; ↗;  
2.5 Ncm; 66 x 16 mm; 63 g

•••► **Hi-Line XS**

consisting parts:  
GB600820 Sterilite spray adaptor Hi

## HiLAN



**GA529** GB559R

### HiLAN with hand control

pneumatic high speed motor for all Hi-Line and Hi-Line XS handpieces  
50 W; 0 - 90.000 1/min; ↗;  
2.5 Ncm; 80 x 19 mm; 132 g

•••► **Hi-Line XS, Hi-Line**

consisting parts:  
GB600820 Sterilite spray adaptor Hi



**GA742R** GB609R

### HiLAN perforator driver

with Hudson chuck  
130 W; 0 - 800 1/min; ↗; 3.6 Nm;  
130 x 37 mm; 495 g  
consisting parts:  
GB600820 Sterilite spray adaptor Hi

## microLAN



**GA553** GB685R

### microLAN 100

pneumatic motor for all micro-Line handpieces  
100 W; 0 - 20.000 1/min; ↗; ↘;  
8.0 Ncm; 95 x 22 mm; 130 g

•••► **micro-Line**

consisting parts:  
GB600810 Sterilite spray adaptor  
micro-Line  
GB600820 Sterilite spray adaptor Hi

# Battery Power Systems

## Acculan 3Ti



**Acculan 3Ti is the battery-operated power system for orthopaedics, traumatology, cardiothoracic surgery, and dermatology**

- > Powerful – up to 250 Watt
- > Reliable – battery and control electronic are not sterilised
- > Easy handling – Titanium housing allows mechanical alkaline reprocessing

Acculan 3Ti – the reinvention of battery power tools



GA671 GB495R

**Acculan 3Ti small drill**

with oscillation and tapping mode;  
not for use with Acetabulum  
reamers  
220 W; ↗; 4.0 mm;  
107 x 47 x 180 mm; 1.100 g  
consisting parts:  
GA666 small NiMH battery  
GA675 lid  
GA678 sterile funnel  
GA672207 K-wire protection sleeve  
GA672211 rinsing adaptor



GA672 GB495R

**Acculan 3Ti drill and reamer**

250 W; ↗; 4.0 mm;  
107 x 47 x 200 mm; 1.250 g  
consisting parts:  
GA675 lid  
GA676 NiMH battery  
GA678 sterile funnel  
GA672207 K-wire protection sleeve  
GA672211 rinsing adaptor



GB635R GB496R

**Acculan 3Ti drill attachment**

with small AO chuck  
0 - 1.000 1/min; 5.0 Nm;  
4.0 mm; 90 x 39 mm; 246 g



GB636R GB496R

**Acculan 3Ti drill attachment**

with Aesculap hexagonal chuck  
0 - 1.000 1/min; 5.0 Nm;  
4.0 mm; 93 x 39 mm; 253 g



GB630R GB496R

**Acculan 3Ti drill attachment**

with Zimmer/Hudson chuck  
0 - 1.000 1/min; 5.0 Nm;  
4.0 mm; 100 x 39 mm; 263 g



GB637R GB496R

**Acculan 3Ti drill attachment**

with trinkle chuck  
0 - 1.000 1/min; 5.0 Nm;  
4.0 mm; 93 x 39 mm; 250 g

Battery

# Acculan 3Ti



**GB638R** **GB496R**  
**Acculan 3Ti drill attachment**  
with Jacobs chuck 0 - 6.5 mm  
0 - 1.000 1/min; 5.0 Nm;  
4.0 mm; 110 x 39 mm; 333 g  
consisting parts:  
GB034800



**GB634R** **GB496R**  
**Acculan 3Ti drill attachment**  
with jacobs chuck keyless 0.6 -  
7.4 mm  
0 - 1.000 1/min; 5.0 Nm;  
4.0 mm; 125 x 39 mm; 430 g



**GB643R** **GB496R**  
**Acculan 3Ti K-wire attachment**  
for K-wires from 0.6 to 4.0 mm  
0 - 1.000 1/min; 93 x 36 x 137 mm;  
300 g



**GB645R** **GB496R**  
**Acculan 3Ti drill attachment**  
for Synthes radiolucent drive  
511.300  
0 - 900 1/min; 110 x 35 mm;  
225 g



**GB654R** **GB496R**  
**Acculan 3Ti reamer attachment**  
with Harris chuck  
0 - 250 1/min; 19.0 Nm;  
4.0 mm; 106 x 39 mm; 323 g



**GB655R** **GB496R**  
**Acculan 3Ti reamer attachment**  
with large AO chuck  
0 - 250 1/min; 19.0 Nm;  
4.0 mm; 107 x 39 mm; 334 g



**GB656R** GB496R  
**Acculan 3Ti reamer attachment**  
with Zimmer/Hudson chuck  
0 - 250 1/min; 19.0 Nm;  
4.0 mm; 118 x 39 mm; 344 g



**GB660R** GB496R  
**Acculan 3Ti sagittal saw attachment**  
for Acculan 3Ti rapid action saw blades with a maximum useable length of 50 mm  
<<>> 0 - 17.000 1/min; 4°47'; 89 x 36 x 44 mm; 270 g



**GA673** GB495R  
**Acculan 3Ti oscillating saw**  
250 W; <<>> 0 - 15.000 1/min;  
 4°47'; 107 x 47 x 200 mm;  
 1.600 g  
consisting parts:  
GA675 lid  
GA676 NiMH battery  
GA678 sterile funnel



**GA668** GB498R  
**Acculan 3Ti straight oscillating saw**  
250 W; <<>> 0 - 15.000 1/min;  
 4°47'; 288 x 54 x 74 mm;  
 1.500 g  
consisting parts:  
GA675 lid  
GA676 NiMH battery  
GA678 sterile funnel



**GA674** GB495R  
**Acculan 3Ti reciprocating saw**  
250 W; <<>> 0 - 15.000 1/min;  
3.2 mm; 180 x 47 x 200 mm;  
 1.540 g  
consisting parts:  
GA675 lid  
GA676 NiMH battery  
GA678 sterile funnel  
GB440215 cleaning plate



**GB436R** GB686R  
**Protection guide**  
for use with GA674, GA874 and saw blade GC761R

# Acculan 3Ti



**GA670**

## Acculan 3Ti dermatome

width keyless adjustable; thickness adjustable (1/10 mm steps); cutting thickness: 0.2 - 1.2 mm; cutting width: max. 78 mm; suitable blades GB228R (to be ordered separately)  
220 W; 0 - 6.500 1/min;  
3.1 mm; 290 x 110 x 90 mm;  
1.700 g  
consisting parts:  
GA666 small NiMH battery  
GA675 lid  
GA678 sterile funnel



**GA677**

## Acculan 3Ti charger

for Acculan 3Ti batteries GA666 and GA676; 4 charging bays; avoids memory effect; automatic battery diagnosis function; with optional available charging bay GA677830 also compatible to Acculan II batteries GA616, GA626 and GA646; power cord not included  
100-240 V AC, 50/60Hz;  
440 x 140 x 180 mm; 3.500 g



**GA666**

## Acculan 3Ti small NiMH battery

only for GA670 and GA671; with integrated electronic motor control  
9,6V ; 1,05Ah  
120 x 43 x 50 mm; 304 g



**GA676**

## Acculan 3Ti NiMH battery

only for GA668, GA672, GA673, GA674; with integrated electronic motor control  
9,6V ; 1,95Ah  
140 x 43 x 50 mm; 425 g



**GA675**

## Acculan 3Ti lid



**GA678**

## Acculan 3Ti sterile funnel



GA679

GB489R

Acculan 3Ti battery removal  
device

Battery



# Handpieces

Hi-Line XS

micro-Line

mini-Line

# Hi-Line XS

## Straight Hi-Line XS handpieces for the Hi-Line XS burr range



GB751R

Hi-Line XS handpiece, straight

1 ring code;  
spray nozzle GB761SU (not  
included)



GB752R

Hi-Line XS handpiece, straight

2 ring code;  
spray nozzle GB762SU (not  
included)



GB753R

Hi-Line XS handpiece, straight

3 ring code;  
spray nozzle GB763SU (not  
included)

## Angled Hi-Line XS handpieces for the Hi-Line XS burr range

Hi-Line XS



**GB756R** GB614R  
Hi-Line XS handpiece, angled  
1 ring code;  
spray nozzle GB761SU (not included)



**GB757R** GB614R  
Hi-Line XS handpiece, angled  
2 ring code;  
spray nozzle GB762SU (not included)



**GB758R** GB614R  
Hi-Line XS handpiece, angled  
3 ring code;  
spray nozzle GB763SU (not included)

# Hi-Line XS

## Angled Hi-Line XS handpieces for the Hi-Line XS burr range



**GB771R** GB614R  
**Hi-Line XS handpiece, angled**  
extra long XL1;  
1 bold ring code;  
spray nozzle GB764SU (not included)



**GB772R** GB614R  
**Hi-Line XS handpiece, angled**  
extra long XL2;  
2 bold ring code;  
spray nozzle GB765SU (not included)



**GB773R** GB614R  
**Hi-Line XS handpiece, angled**  
extra long XL3;  
3 bold ring code;  
spray nozzle GB766SU (not included)

## Craniotome and multipurpose handpiece



**GB740R** GB614R  
**Hi-Line XS craniotome and multi purpose handpiece**  
 to be used with dura guards, drill depth guard and holding sleeve  
 GB741R - GB748R,  
 spray nozzle GB767SU (not included)



**GB742R** GB615R  
**Hi-Line XS dura guard**  
 fixed;  
 2 ring code;  
 to be used with craniotome burr  
 GE520R, GE520SU, GE529SU



**GB746R** GB615R  
**Hi-Line XS dura guard**  
 steerable;  
 1 ring code;  
 to be used with craniotome burr  
 GE420R, GE420SU, GE429SU



**GB741R** GB615R  
**Hi-Line XS dura guard**  
 fixed;  
 1 ring code;  
 to be used with craniotome burr  
 GE420R, GE420SU, GE429SU



**GB743R** GB615R  
**Hi-Line XS dura guard**  
 fixed;  
 3 ring code;  
 to be used with craniotome burr  
 GE620R, GE620SU, GE629SU



**GB747R** GB615R  
**Hi-Line XS dura guard**  
 steerable;  
 2 ring code;  
 to be used with craniotome burr  
 GE520R, GE520SU, GE529SU

# Hi-Line XS



**GB748R** GB615R

**Hi-Line XS dura guard**  
steerable;  
3 ring code;  
to be used with craniotome burr  
GE620R, GE620SU, GE629SU



**GB744R** GB615R

**Hi-Line XS drill depth guard**  
for craniotome handpiece and  
multipurpose handpiece GB740R  
2 - 8 mm



**GB745R** GB615R

**Hi-Line XS holding sleeve**  
for craniotome handpiece and  
multipurpose handpiece GB740R



# micro-Line

For burrs with shank diameter of 2.35 mm



**GB109** **GB681R**

**micro-Line handpiece**

coding: 1 blue ring;

for 51 mm burrs

1:1; 81 mm



**GD450M** **GB681R**

**micro-Line handpiece**

coding: 1 blue ring;

for 70 mm burrs;

suitable spray nozzle GD460R (not included)

1:1; 97 mm



**GD451M** **GB681R**

**micro-Line handpiece**

coding: 1 blue ring;

for 95 mm burrs;

1:1; 127 mm

For burrs with shank diameter of 2,35 mm; angulation 21°



**GD455M** GB681R  
**micro-Line handpiece, angled**  
coding: 1 red ring;  
for 70 mm burrs;  
suitable spray nozzle GD460R (not  
included)  
1:2; 120 mm



**GD456M** GB681R  
**micro-Line handpiece, angled**  
coding: 1 red ring;  
for 95 mm burrs;  
1:2; 150 mm



**GD457M** GB681R  
**micro-Line handpiece, angled**  
coding: 1 red ring;  
for 125 mm burrs;  
1:2; 180 mm

micro-Line

# micro-Line

For use with micro-inserts with imperator gearing, angulation 21°



GB114R GB681R

**micro-Line handpiece**

coding: 1 red ring;  
suitable spray nozzle GA227R (not  
included)

1:2; 77 mm



GB111 GB681R

**micro-Line handpiece**

coding: 1 green ring;  
suitable spray nozzle GA227R (not  
included)

4:1; 74 mm



GB113 GB681R

**micro-Line handpiece**

coding: 2 green rings;  
suitable spray nozzle GA227R (not  
included)

30:1; 74 mm



GB173 GB690R

**Imperator attachment**

for burrs with 2.35 mm diameter; up  
to 70 mm length; minimal shaft  
length 30 mm; spray nozzle GA227R  
(not included)

1:1; 43 mm



GB164

**Mounting tool**

for imperator attachment GB173



GB119

**Contra angle head**

coding: 1 green dot;  
for tools with Dental shank

2:1; 43 mm



**GD305** **GB681R**  
**micro-Line reciprocating saw**  
keyless; suitable spray nozzle  
GD306R (not included)  
2,5:1;  $\lll \ggg$  8.000 1/min;  $\downarrow$  1.3 mm;  
 $\blacksquare$  48 g  
consisting parts:  
GD320R saw blade



**GD307** **GB681R**  
**micro-Line sagittal saw**  
with captive nut; suitable spray  
nozzle GD308R (not included)  
1:1;  $\lll \ggg$  20.000 1/min;  $\swarrow$  7°;  
 $\blacksquare$  50 g  
consisting parts:  
GD372206 key



**GB129** **GB682R**  
**micro-Line sagittal saw**  
spray nozzle GA222R (not included)  
1:1;  $\lll \ggg$  20.000 1/min;  $\swarrow$  4°46';  
 $\blacksquare$  345 g  
consisting parts:  
TE472 ring spanner  
GB125804 nut



**GB128R** **GB682R**  
**micro-Line transversal saw**  
suitable spray nozzle GA221R (not  
included)  
1:1;  $\lll \ggg$  20.000 1/min;  $\swarrow$  10°24';  
 $\blacksquare$  385 g  
consisting parts:  
GB020R allen key



**GB130R** **GB682R**  
**micro-Line reciprocating saw**  
suitable spray nozzle GA223R (not  
included)  
1:1;  $\lll \ggg$  20.000 1/min;  $\downarrow$  3.0 mm;  
 $\blacksquare$  310 g



**GB200** **GB682R**  
**micro-Line drill pistol**  
spray nozzle GA221R (not included)  
 $\blacktriangleleft \blacktriangleright$ ; 20:1;  $\blacksquare$  2.7 mm;  $\blacksquare$  700 g

# micro-Line



**GB222**

**Jacobs chuck**

Ø 0 - 4 mm

consisting parts:  
GA062R Key



**GB223**

**Jacobs chuck**

Ø 0 - 6.5 mm

consisting parts:  
GA031R key



**GB224**

**Small AO chuck**



**GB227R**

**K-wire driver**

for K-wires from 0.6 to 2.5 mm



**GB231R**

**GB682R**

**micro-Line Dermatome acc. to  
Wagner**

cutting width keyless adjustable;  
thickness adjustable (1/10 mm  
steps);  
cutting thickness: 0.2 - 1.2 mm;  
cutting width: max. 78 mm  
6:1; 3.0 mm; 770 g  
consisting parts:  
GB228R blades



**GB169R**

**GB682R**

**micro-Line drill**

for exchangeable chucks; suitable  
spray nozzle GA228R (not included)  
20:1; 280 g  
consisting parts:  
GC523200 ring spanner



**GB106R**  
**Hudson chuck**  
for GB169R

**GB268R**    **GB682R**  
**micro-Line craniotome**  
with fixed duraguard  
1:1;  $\ddot{m}$  200 g  
consisting parts:  
GB267R handpiece  
GB298R fixed dura-guard  
GC303R burr



**GB265R**    **GB682R**  
**micro-Line craniotome**  
with turnable duraguard  
1:1;  $\ddot{m}$  200 g  
consisting parts:  
GB266M turnable dura-guard  
GB267R handpiece  
GC303R burr



**GB292**    **GB681R**  
**micro-Line pediatric craniotome,**  
**Freiburg pattern**  
with fixed duraguard  
1:1;  $\ddot{m}$  100 g  
consisting parts:  
GB109 handpiece  
GB291R dura-guard  
GC304R burr

# mini-Line

## mini-Line



**GB376R** **GB673R**  
**mini-Line burr attachment**  
straight; for 2,35 mm shank; for  
burrs up to 70 mm length; suitable  
spray nozzle GB378 (not included)  
1:2; 108 g



**GB382R** **GB673R**  
**mini-Line drill attachment**  
for tools with AO shank; suitable  
spray nozzle GB378 (not included)  
16:1; 98 g



**GB383R** **GB673R**  
**mini-Line drill attachment**  
Jacobs chuck; suitable spray nozzle  
GB378 (not included)  
Ø 0 - 4 mm  
16:1; 2.5 mm; 160 g  
consisting parts:  
GA062R Key



**GB377R** **GB674R**  
**mini-Line heavy duty drill  
attachment**  
Jacobs chuck; suitable spray nozzle  
GB378 (not included)  
Ø 0 - 6.5 mm  
36:1; 2.5 mm; 220 g  
consisting parts:  
GA031R key



**GB389R** **GB674R**  
**mini-Line drill attachment**  
with quick coupling for tools  
with hexagonal shank; suitable  
spray nozzle GB378 (not included)  
36:1; 2.5 mm; 145 g



**GB373R** **GB673R**  
**mini-Line angled drill attachment**  
for tools with dental shank  
8:1; 91 g



**GB398R** **GB674R**  
**mini-Line K-wire attachment**  
for use with GD684;  
for K-wires from 0.6 to 2.5 mm  
16:1; 190 g



**GB390R** **GB673R**  
**mini-Line sagittal saw attachment**  
1:1; 5°; 100 g



**GB392R** **GB674R**  
**mini-Line transversal saw attachment**  
1:1; 10°; 140 g



**GB374R** **GB673R**  
**mini-Line angled K-wire attachment**  
for use with straight mini motors;  
for K-wires from 0.6 to 1.6 mm  
17:1; 120 g



**GB391R** **GB673R**  
**mini-Line reciprocating saw attachment**  
1:1; 2.1 mm; 100 g



**GA187R** **GB673R**  
**micro-Line adaptor**  
handpiece: micro;  
motor: mini

# Accessories

## Accessories for electric motors



**GA411**

**microspeed uni mobile stand**  
adjustable height  
1200 - 1500 mm



**GA415**

**Mobile stand**  
for Elan, microspeed EC and  
Microtron devices; adjustable height  
1200 - 1500 mm



**GD416M**

**Mobile stand**  
for AESCULAP control units  
1320 mm  
consisting parts:  
GD414R



**GD412804**

**Bottle holder**



**GA414**

**Sterile disposable drape**  
for GD417, GA419 and GD421  
②; package of 50 pieces



**TE780**  
**Power cord Europe; black**

1.5 m

**TE730**  
**Power cord Europe; black**

5.0 m

**TE734**  
**Power cord UK; black**

5.0 m

**TE735**  
**Power cord USA, Canada, Japan;  
black**

3.5 m

# Accessories


Wall outlet	Motor inlet	Air hoses
		<b>GA505R</b>  <b>Air hose</b> 5.0 m
AESCULAP Dräger	small	<b>GA513R</b>  <b>Air hose</b> 3.0 m
Wall outlet	Motor inlet	
		<b>GA523R</b>  <b>Air hose</b> (for irrigation pump GA530) 2.3 m
AESCULAP Dräger	small	
Wall outlet	Motor inlet	
		<b>GA507R</b>  <b>Air hose</b> 5.0 m
DIN	small	

Wall outlet	Motor inlet	
		<b>GA506R</b> <b>GB575R</b> Air hose  5.0 m
Schrader with diffusor	small	
		<b>GA464R</b> <b>GB463R</b> <b>GA466R</b> Air hose  3.0 m  5.0 m
AESCULAP Dräger	large	
		<b>GA460R</b> <b>GB463R</b> Air hose, spiral type  3.0 m

Recommended pressure for all pneumatic AESCULAP Systems is 6 - 10 bar, 85 - 145 psi  
 Motor connection small fits: GA502, GA508R, GA519, GA520R, GA529, GA551, GA552, GA553, GA554, GA740R, GA741R, GA742R  
 Motor connection large fits: GA200, GA430, GA437, GA440, GA521, GA530

# Accessories

Wall outlet	Motor inlet		
		<b>GA463R</b> Air hose  3.0 m	<b>GA465R</b> Air hose  5.0 m
DIN	large		
		<b>GA461R</b> Air hose  3.0 m	
DIN with diffusor	large		
		<b>GA468R</b> Air hose  5.0 m	<b>GB480R</b>
Schrader with diffusor	large		
<p>Recommended pressure for all pneumatic AESCULAP Systems is 6 - 10 bar, 85 - 145 psi Motor connection small fits: GA502, GA508R, GA519, GA520R, GA529, GA551, GA552, GA553, GA554, GA740R, GA741R, GA742R Motor connection large fits: GA200, GA430, GA437, GA440, GA521, GA530</p>			

Wall outlet	Motor inlet	Extension hoses
		<b>GA194R</b> Extension hose  3.0 m
AESCLAP Dräger	AESCLAP Dräger	<b>GA195R</b> Extension hose  5.0 m
Wall outlet	Motor inlet	
		<b>GA087</b> Extension hose  5.0 m
DIN	DIN	
Pneumatic accessories		
	<b>GA521</b> <b>Foot control</b> for pneumatic systems air inlet: large air outlet: AESCLAP Dräger	<b>GA099</b> <b>Air pressure reducer</b> air pressure reducer with DIN quick lock coupling; for compressed air cylinders acc. to DIN477 G5/8" AG; air inlet 100 - 200 bar / 2130-2840 psi consisting parts: GA088C DIN quick lock coupler

# Accessories



**GA088C**

**Quick lock coupler**  
for air pressure reducer GA099; for  
DIN system



**GA508R**

**Angled hose adaptor**  
for all small Aesculap motors



**GA208**

**GB474R**

**Adaptor for pneumatic systems**  
Synthes motor - air hose with  
Aesculap motor connector large



**GA091C**

**Quick lock coupler**  
for french system



**GA206R**

**Adaptor for pneumatic systems**  
AESCULAP motor large - Synthes air  
hoses

## Irrigation

**GA530****Irrigation pump**

for pneumatic systems; to be used with tubing set GA391SU and GA526P (not included)  
pump flow 0– 80 ml/min  
100–240 V AC, 50/60Hz;  
305 x 305 x 135 mm; 5.500 g  
consisting parts:  
GD412804 bottle holder

**GA392SU****Tubing set with Y-Adaptor**

complete;  
including adaptor for spray nozzles,  
fixations clips and punch for  
irrigation bottle;  
for microspeed uni  
⑧; ⑩ package of 10 pieces

**GA526P****Fixation clips**

for combining tubing set GA391 or  
GA391SU with pneumatic motors  
⑩ package of 8 pieces

**GA391SU****Tubing Set**

complete;  
including adapter for spray nozzles,  
fixation clips and punch for  
irrigation bottle;  
for microspeed uni, Microspeed EC,  
Microtron EC, HiLAN (in connection  
with GA526P)  
⑧; ⑩ package of 10 pieces

**GA391****Tubing set**

complete; including adaptor for  
spray nozzles, fixation clips and  
punch for irrigation bottle; for Elan  
(micro and mini cable only); HiLAN  
in connection with GA526P  
⑧; ⑩ package of 10 pieces

# Accessories

## Adaptor



GB421R    GB467R

### Hudson adaptor

tool: Hudson;  
motor: Harris

3.8 mm



GB422R    GB467R

### AO large adaptor

tool: AO large;  
motor: Harris

3.8 mm



GB423R    GB467R

### DIN adaptor

tool: DIN;  
motor: Harris

3.8 mm



GB425R    GB467R

### Zimmer/Hudson adaptor

tool: Zimmer/Hudson;  
motor: Harris

3.8 mm



GB426R    GB467R

### Jacobs chuck

tool: Jacobs;  
motor: Harris

Ø 0 - 6,5 mm

3.8 mm



GB628R

### Synthes DHS/DCS adaptor

tool: Synthes DHS/DCS;  
motor: AO large

4.2 mm

GB182R GB690R**Harris adaptor**

tool: Harris;  
motor: Hudson  
not canulated

GB184R GB690R**Harris adaptor**

tool: Harris;  
motor: AO large  
 3.5 mm

GB178R GB690R**Harris adaptor**

tool: Harris;  
motor: DIN  
 3.5 mm

XF457R GB690R**DIN adaptor**

tool: DIN;  
motor: AO large  
 3.8 mm

GB007R GB690R**AO large adaptor**

tool: AO large;  
motor: DIN  
 4.0 mm

GB187R GB690R**Harris adaptor**

tool: Harris;  
motor: triangle  
 3.5 mm

# Accessories



GA187R [GB673R]

**micro-Line adaptor**

handpiece: micro;  
motor: mini

## Keys



GB034R [GB690R]

**Key**

suitable for GB223R, GB386R,  
GB416R, GB426R, GB424R, GB177R,  
GB638R



GA186209 [GB690R]

**Key**

for mini-coupling  
suitable for GA186, GA551, GA552,  
GA647, GD674, GD684



TE472 [GB690R]

**Ring spanner**

for Aesculap star coupling



GB020R [GB690R]

**Allen Key**

suitable for GB101, GB102R,  
GB115R, GB121R, GB125R, GB210R,  
GB230R, GB377R



**GA062R**    **GB690R**

**Key**

for small 4mm Jacobs chuck  
suitable for GB108R, GB166R,  
GB181R, GB222R, GB383R, GB415R



**GD372206**

**Key**

width across flats 6,0 mm  
suitable for GD307



**GC523R**

**Ringspanner**

suitable for GA522R, GB106R,  
GB108R, GB169R

## Skin graft mesher



**BA720R**    **GB688R**

**Skin graft mesher**

for use in combination with skin carriers BA721 (1:1,5) / BA722 (1:3) / BA723 (1:6)

150 x 195 x 105 mm; 4.600 g  
consisting parts:  
BA725R cutting roller  
BA726R ratchet wrench



**BA725R**

**Cutting roller**

# Accessories

## Cement injection



**GB141**

### Cement syringe

for cement injector GB171N; with special injection tubing

50 ml

②: 10 package of 10 pieces



**GB171N**

### Cement injector w. Harris shank

for disposable syringe GB141

## Bone mill



**GB040**

### Bone mill

for use with motor or crank handle consisting parts:

GB042R crushing drum



**GB042R**

### Crushing drum

standard; Harris shaft; section of the produced bone chip 0.7 x 3 mm



**GB043R**

### Crushing drum

fine; Harris shaft; section of the produced bone chip 0.7 x 2 mm



**GB044R**

### Crushing drum

coars; Harris shaft; section of the produced bone chip 0.7 x 4 mm



**GB041R**  
Crank handle for GB040



**GB046**  
Insert for GB040  
for small bone grafts



# Processing

## Sterilit Power Systems

### Eccos

- Eccos for Power Systems
- Eccos for handpieces
- Eccos for accessories
- Eccos sets

# Sterilit Power Systems

## Sterilit Power Systems



**GB600**  
**Sterilit Power Systems**  
300 ml  
6 pieces



**GB600820**  
**Sterilit Power Systems spray adaptor**  
for HiLAN, HiLAN XS, Hi-Line, Hi-Line XS, Microspeed EC, microspeed uni, microLAN, miniLAN



**GB600840**  
**Sterilit Power Systems spray adaptor**  
for Acculan 3Ti



**GB600810**  
**Sterilit Power Systems spray adaptor**  
for micro-Line



**GB600830**  
**Sterilit Power Systems spray adaptor**  
for mini-Line



**GB600870**  
**Sterilit Power Systems spray adaptor**  
for microspeed uni XS



**GA059**

**Sterilit drip lubricator**

50 ml

# Eccos

## Eccos

### **The solution for storage and mechanical reprocessing of your Aesculap Power Systems**

Basic requirements for mechanical reprocessing in combination with the Eccos system:

- > A suitable cleaner must be used (pH < 11 for Aluminum components)
- > The temperature in the cleaning cycle must not exceed 60°C
- > Clean for at least 3 minutes
- > Apply intermediate rinse for at least 1 minute
- > Carry out intensive final rinse with distilled, demineralised or fully desalinated water
- > Only Thermal disinfection: 5 minutes rinse at max. 96°C with distilled, demineralised or fully desalinated water
- > Complete the programme with a drying phase of at least 10 minutes at a temperature of not more than 120°

## microspeed uni/microspeed arthro



**GB579R**

Eccos fixation

for GD684

56 mm



**GB689R**

Eccos fixation

for GD672, GD673

56 mm



**GB587R**

Eccos fixation

for GD685

56 mm



**GB578R**

Eccos fixation

for GD674, GD676, GD678

56 mm



**GB580R**

Eccos fixation

for GD682, GD681

56 mm



**GB586R**

Eccos fixation

for GD675

56 mm



# Eccos for Power Systems

## Elan EC



**GB680R**

**Eccos fixation**  
for GA173, GA176

56 mm



## HiLAN XS/HiLAN



**GB577R**

**Eccos fixation**  
for small Aesculap motors with  
mounted angled hose adaptor  
GA508R

56 mm



**GB610R**

**Eccos fixation**  
for GA740R, GA741R

56 mm



**GB609R**

**Eccos fixation**  
for GA742R, GA522R

56 mm



## microLAN



**GB685R**

### Eccos fixation

for GA502, GA553, GA554, GD657,  
GD658

56 mm



## Acculan 3Ti



**GB495R**

### Eccos fixation

for GA671, GA672, GA673, GA674,  
GA871, GA873, GA874

106 mm



**GB496R**

### Eccos fixation

for GB630R, GB634R, GB635R,  
GB636R, GB637R, GB638R, GB643R,  
GB644R, GB645R, GB654R, GB655R,  
GB656R, GB660R

106 mm



**GB497R**

### Eccos fixation

for GB630R, GB634R, GB635R,  
GB636R, GB637R, GB638R, GB643R,  
GB644R, GB645R, GB654R, GB655R,  
GB656R, GB660R

106 mm



**GB487R**

### Eccos fixation

for GA675

56 mm



# Eccos for Power Systems



**GB488R**  
**Eccos fixation**  
for GA678  
76 mm  
 76 mm



**GB489R**  
**Eccos fixation**  
for GA614, GA644, GA679  
106 mm  
 106 mm



**GB498R**  
**Eccos fixation**  
for GA668, GA670, GA870  
106 mm  
 106 mm

## Acculan 3Ti Electro



**GB485R**  
**Eccos fixation**  
for GA875  
56 mm  
 56 mm



**GB486R**  
**Eccos fixation**  
for GA878  
56 mm  
 56 mm

## Hi-Line XS



**GB611R**

**Eccos fixation for 1  
Hi-Line XS handpiece**

for GB740R, GB751R, GB752R,  
GB753R, GB756R, GB757R, GB758R,  
GB771R, GB772R, GB773R

56 mm




**GB613R**

**Eccos fixation**

for 13 burs of the Hi-Line XS, Hi-Line or 2,35 mm shaft range

36 mm



**GB614R**

**Eccos fixation for 3  
Hi-Line XS handpieces**

for GB740R, GB751R, GB752R,  
GB753R, GB756R, GB757R, GB758R,  
GB771R, GB772R, GB773R

56 mm



**GB615R**

**Eccos fixation**

for GB741R - GB748R

56 mm



## micro-Line



**GB681R**

**Eccos fixation**

for GB109, GB110, GB111, GB112,  
GB113, GB114R, GB291R, GB292R,  
GD305, GD307, GB310, GD450M,  
GD451M, GD452M, GD455M,  
GD456M, GD457M, GD465M,  
GD466M, GD467M

56 mm




**GB686R**

**Eccos fixation**

for GB266, GB298, GB436R

56 mm



# Eccos for handpieces



**GB682R**

**Eccos fixation**

for GB128R, GB129, GB130R,  
GB166R, GB169R, GB200, GB231R,  
GB265R, GB267R, GB268R, GB280,  
GB801R

76 mm



## mini-Line



**GB673R**

**Eccos fixation**

for GB373R, GB374R, GB376R,  
GB377R, GB379R, GB380R, GB381R,  
GB382R, GB383R, GB386R, GB390R,  
GB391R

56 mm



**GB674R**

**Eccos fixation**

for GB385R, GB388R, GB389R,  
GB392R, GB397R, GB398R

56 mm



## Skin graft mesher



**GB688R**  
**Eccos fixation**  
for BA720R  
106 mm

## Accessories



**GB564R**  
**Eccos fixation**  
for GA505R, GA513R, GA523R  
56 mm



**GB463R**  
**Eccos fixation**  
for GA460R, GA464R, GA466R  
76 mm



**GB575R**  
**Eccos fixation**  
for GA506R  
76 mm



**GB576R**  
**Eccos fixation**  
for GA507R  
56 mm

# Eccos for accessories



**GB480R**  
**Eccos fixation**  
for GA468R  
76 mm



**GB474R**  
**Eccos fixation**  
for GA206R  
56 mm




**GB572R**  
**Eccos fixation**  
for GB300R - GB308R  
56 mm




**GB467R**  
**Eccos fixation**  
for GB421R, GB422R, GB423R,  
GB424R, GB425R  
56 mm



**GB481R**  
**Eccos fixation**  
for saw blades and other accessories  
56 mm



**GB690R**  
**Eccos fixation**  
for GB173R, GB119, GB154, GB194,  
GB628R, GA031201, GB424201,  
GA437224, GB009C  
36 mm





**GB459R**  
**Eccos Montageset**

# Eccos sets


## Acculan 3Ti



**GB243R**

**Eccos Set for 2 Acculan 3Ti machines**

106 mm  
consisting parts:  
GB243800  
GB481R  
GB487R  
GB488R  
GB495R  
GB496R  
GB497R  
GB690R



**GB244R**

**Eccos Set for 1 Acculan 3Ti machine**

106 mm  
consisting parts:  
GB243800  
GB481R  
GB487R  
GB488R  
GB495R  
GB496R  
GB497R  
GB690R



**GB256R**

**Eccos Set for Acculan 3Ti dermatome**

106 mm  
consisting parts:  
GB487R  
GB488R  
GB498R  
JC254R

## microspeed uni




**GB255R**

**Eccos Set for microspeed uni**

Neurosurgery

56 mm  
consisting parts:  
GB586R  
GB587R  
GB614R  
GB615R  
GB689R  
JF212R



**GB254R**

**Eccos Set for microspeed uni**

small bone surgery

56 mm  
consisting parts:  
GB579R  
GB673R  
GB674R  
GB690R  
JF212R

## HiLAN XS



**GB252R**

**Eccos Set for HiLAN XS**

56 mm  
consisting parts:  
GB564R  
GB609R  
GB610R  
GB614R  
GB615R  
JF213R



# Index

# Index by reference numbers

BA720R	57	GA676	28	GB243R	74	GB587R	65
BA725R	57	GA677	28	GB244R	74	GB600	62
		GA678	18,28	GB252R	75	GB600810	62
GA059	63	GA679	29	GB254R	75	GB600820	62
GA062R	57	GA740R	23	GB255R	75	GB600830	62
GA087	51	GA742R	23	GB256R	74	GB600840	62
GA088C	52	GA830	21	GB265R	43	GB600870	62
GA091C	52	GA835	21	GB268R	43	GB609R	66
GA099	51	GA870	17	GB292	43	GB610R	66
GA173	21	GA871	17	GB373R	44	GB611R	69
GA176	21	GA873	17	GB374R	45	GB613R	69
GA186209	56	GA874	17	GB376R	44	GB614R	69
GA187R	45,56	GA875	18	GB377R	44	GB615R	69
GA188	21	GA876	18	GB382R	44	GB628R	54
GA194R	51	GA877	17	GB383R	44	GB630R	25
GA195R	51	GA878	18	GB389R	44	GB634R	26
GA206R	52			GB390R	45	GB635R	25
GA208	52	GB007R	55	GB391R	45	GB636R	25
GA391	53	GB020R	56	GB392R	45	GB637R	25
GA391SU	53	GB034R	56	GB398R	45	GB638R	26
GA392SU	53	GB040	58	GB421R	54	GB643R	26
GA411	46	GB041R	59	GB422R	54	GB645R	26
GA414	46	GB042R	58	GB423R	54	GB654R	26
GA415	46	GB043R	58	GB425R	54	GB655R	26
GA460R	49	GB044R	58	GB426R	54	GB656R	27
GA461R	50	GB046	59	GB436R	17,27	GB660R	27
GA463R	50	GB106R	43	GB459R	73	GB673R	70
GA464R	49	GB109	38	GB463R	71	GB674R	70
GA465R	50	GB111	40	GB467R	72	GB680R	66
GA466R	49	GB113	40	GB474R	72	GB681R	69
GA468R	50	GB114R	40	GB480R	72	GB682R	70
GA505R	48	GB119	40	GB481R	72	GB685R	67
GA506R	49	GB128R	41	GB485R	68	GB686R	69
GA507R	48	GB129	41	GB486R	68	GB688R	71
GA508R	52	GB130R	41	GB487R	67	GB689R	65
GA513R	48	GB141	58	GB488R	68	GB690R	72
GA521	51	GB164	40	GB489R	68	GB740R	35
GA523R	48	GB169R	42	GB495R	67	GB741R	35
GA526P	53	GB171N	58	GB496R	67	GB742R	35
GA529	23	GB173	40	GB497R	67	GB743R	35
GA530	53	GB178R	55	GB498R	68	GB744R	36
GA553	23	GB182R	55	GB564R	71	GB745R	36
GA666	28	GB184R	55	GB572R	72	GB746R	35
GA668	27	GB187R	55	GB575R	71	GB747R	35
GA670	28	GB200	41	GB576R	71	GB748R	36
GA671	25	GB222	42	GB577R	66	GB751R	32
GA672	25	GB223	42	GB578R	65	GB752R	32
GA673	27	GB224	42	GB579R	65	GB753R	32
GA674	27	GB227R	42	GB580R	65	GB756R	33
GA675	28	GB231R	42	GB586R	65	GB757R	33

# Index by reference numbers

GB758R	33
GB771R	34
GB772R	34
GB773R	34
GC523R	57
GD305	41
GD307	41
GD372206	57
GD412804	46
GD416M	46
GD450M	38
GD451M	38
GD455M	39
GD456M	39
GD457M	39
GD665	13
GD668	9
GD670	9
GD671	9,13
GD672	9,13
GD673	9,13
GD674	10,14
GD675	10
GD676	9
GD678	10,14
GD681	10,13
GD682	11,13
GD684	10,14
GD685	10
TE472	56
TE730	47
TE734	47
TE735	47
TE780	47
XF457R	55



The main product mark 'Aesculap' is a registered mark of Aesculap AG.

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 (0) 74 61 95-0 | Fax +49 (0) 74 61 95-26 00 | [www.aesculap.com](http://www.aesculap.com)

Aesculap – a B. Braun company

Subject to technical changes. All rights reserved.  
This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.

Brochure No. 022702

0111/3.5/3